PROTOCOL

"Anchored Cell detachment, concentration and cell wash for further studies"

Without Celartia Centrifuge

PETAKA G3 DESCRIPTION

	Petaka G3		
Biotechnology	New. Patented		
	Self regulated gas diffusion system.		
Air space (volume lost)	3.8 % of the total volume		
Surface available for cell growth	150 cm2		
Average Cell Production	20 million cells		
Device volume	59 cm3		
Design	 New. Patented. Adapted to robotic cell culture Prepared for internal cell concentration by centrifugation. Prepared to segregate cells from supernatant. Prepared for direct positive and negative immune-magnetic cell sorting Prepared for Microgravity cell culture (orbital incubation) Prepared for long distance cell shipping without freezing or cooling. Prepared for induction of "in vitro cell dormancy" 		

PROTOCOL. Adherent Cell Transfer without Petaka Centrifuge

Materials needed

PetakaG3 Flame 30 mL Syringe Petaka tips 10 mL conical tubes Laminar flow hood Centrifuge Cell counter

27% savings in instrumentation

Disinfect Petaka port with 90 &% Ethanol and flame it briefly Time:0.1 min

Completely withdraw the media Time:1 min

Inject 4 mL of Trypsin solution Time:0.1 min

Rotate and gently shake the PETAKA to cover all cell culture surface with a thin layer of trypsin Time: 1 min

Incubate for 3 to 4 minutes at 37°C Time 4 min



Check cell detachment under the microscope Time: 1 min	
Disinfect Petaka port with 90% Ethanol and flame it briefly Time:0.1 min	
Inject into Petaka 6 mL of media having 20% FBS to inhibit the trypsin activity Time: 0.1 min	
When the injection is complete, withdraw a small sample of cell suspension Time: 0.5min	
Count the cell content in cells/mL and estimate the total amount of cells during the first centrifugation Time:5 min compensated by centrifugation time	
Completely withdraw the media Time:1 min	
IN LAMINAR FLOW HOOD. Transfer cell suspension to a 10 mL conical tube Close the tube with cap. Time: 1 min	
Close the tube with the cap	
Centrifuge the conical tube at the necessary RPM and for sufficient time to obtain a soft cell pellet Time: 6 mins	
IN LAMINAR FLOW HOOD. Decant the media keeping intact the cell pellet. Time: 0.1 min	
IN LAMINAR FLOW HOOD. Add 10 mL of PBS to the conical tube and vortex it until complete cell suspension. Time: 2 mins	
Close the tube with the cap	

Centrifuge the conical tube at the necessary RPM and time to obtain a soft cell pellet Time: 6 mins	
IN LAMINAR FLOW HOOD. Decant the media keeping the cell pellet intact . Time: 0.1 min	
Transfer the cells to the desired device adjusting volumes	-4
TOTAL Time: 24.1mins 15.7% savings time compared to T75 flask	